


QUANTITATIVE EVALUATION WITH *visionCATS*

visionCATS controls all functions of the TLC Scanner 4 and enables quantitative evaluation of the generated densitometric data. Only positioning of the object to be measured is performed manually and, if desired, switching on the internal illumination to assist correct positioning. Optimal settings of the electronic amplification are automatically selected for scanning in absorbance and/or fluorescence mode respectively. The 16 bit A/D converter ensures optimally adapted resolution of the measurement signal.


To determine the substance concentration in a sample, five different quantification functions (e.g. linear and polynomial) are available. Several scanning steps (e.g. scanning the plate after development and scanning the same plate after derivatization) and up to five different evaluations can be performed (with data obtained from single wavelength, multiple wavelengths or a combination of measurements in absorbance and fluorescence detection mode).

The Scanner Ultimate Package:


- **Multi Wavelength Scan:** this feature offers the possibility to perform a multi-wavelength scan with up to 31 selected wavelengths or a combination of measurements in absorbance and fluorescence detection mode.
- **Scanner Quantification:** this feature allows to quantify each individual substance on the plate. Five different quantification functions are available for evaluation to determine the concentration of the substance in a sample. In one analysis file up to five evaluation steps can be performed in multiple plate states (e.g. plate after development and same plate after derivatization).
- **Spectrum Scan:** this feature includes the measurement of the spectrum of each individual substance on the plate including the evaluation of the substance purity by comparison with reference standard.



Densitograms are displayed in 3D, top or front view. Several peak integration and baseline correction settings can be selected.



The separated compounds are assigned. For quantification, data from the multi-wavelength scan at the optimum wavelength for each compound is used.



For evaluation the best fitting calibration model is used. Quantitation can be done via peak height or area.

Ordering Information

027.6200 CAMAG® TLC Scanner 4

for the densitometric evaluation of TLC/HPTLC plates, spectral range 190 to 900 nm, plate sizes up to 20 x 20 cm, absorbance and fluorescence mode, CAMAG *visionCATS* controlled.

028.0000 CAMAG® HPTLC Software *visionCATS* Basic Version

including access and control of all instruments - one server, one client, Instrument Diagnostics (xQ), analytical reports, access to Method Library. Needs to be purchased separately and is not included in any Ultimate Package.

028.3000 CAMAG® HPTLC Software *visionCATS*: Scanner Ultimate Package

Multi Wavelength Scan, Spectrum Scan and Quantification.

CAMAG HPTLC Software *visionCATS* Basic Version (028.0000) needs to be purchased separately.

Before installing *visionCATS*, please visit www.camag.com/visionCATS for recommended system requirements and further information.