Ultra Yield™ Flasks – Brief Instructions

Culture medium

Poor media like LB (Luria Broth) cannot generate more biomass. **Use enriched media like TB (Terrific Broth) or modified TB like SRAM:** yeast extract (24 g), tryptone (16 g), casamino acids (10 g), 1 % glycerol w/v and 100 mM MOPS buffer (pH 7.4) per L of media. Adapt buffer pH to the bacterial strain of interest. Incubate for 30 h, up to 7 days. As a guideline for *E. coli*, an OD₆₀₀ of up to 30 can be expected (depending on the clone).

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Shaking speed

Ultra Yield™ Flasks can provide 10x increased aeration over standard shake flasks, but usual shaking speeds of 200–220 min⁻¹ (25 mm shaking throw) are not enough. **Use minimum 300 min⁻¹, preferably up to 400 min⁻¹ (25 mm shaking throw). 2.5 L flask with 1 L filling volume: 250 min⁻¹.** INFORS HT recommends to use Ultra Yield™ Flasks only on 25 mm shaking diameter. 50 mm: reduce min⁻¹ by 40 %.

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Filling volume

Ultra Yield™ Flasks can be used with up to 40 % of the nominal volume. Maximum gas transfer is achieved with about 20 %. Perform the first experiments with the optimum filling volume to see how much biomass you can get. Then, use the maximum filling volume, which may lead to less biomass. **Note that choosing the optimum volume may greatly reduce formation of secondary proteins, saving costs and time in your subsequent protein purification steps.**

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