

VIP

Горизонтальный морозильник-80°C

84 л



MDF-C8V1-PE



Горизонтальный морозильник PHCbi MDF-C8V1 VIP -80 ° C обеспечивает максимальную емкость для хранения образцов при оптимальной занимаемой площади в сочетании с лучшими в отрасли показателями охлаждения и надежности для идеального решения долгосрочного хранения.

Оптимальное хранение образцов


Использование запатентованных панелей VIP PLUS в морозильнике обеспечивает уменьшенную толщину стенок для максимального внутреннего объема. Эта технология позволяет MDF-C8V1 занимать меньше места для установки при большой емкости.

Надежные технологии PHCbi

Компрессор Cool Safe®, специально разработанный для применения при сверхнизких температурах, используется в проверенной PHCbi автоматической каскадной системе охлаждения, обеспечивающей высочайший уровень производительности и надежности.

Простота использования и интеллектуальная безопасность


Комплексные функции установки, сигнализации, мониторинга и диагностики, контролируемые встроенным микропроцессорным контроллером PHCbi с цифровым отображением всех функций.



Personal Storage Ideal for reliable storage of small volume clinical trial samples and temperature control testing.



Dependable Preservation Quality of design and manufacture ensures trusted, reliable storage that maintains the integrity of precious samples.



Uniform Sample Storage The uniform and stable temperatures provided are ideal for storage of DNA, RNA and other important biological samples.

