


Same Foot Print Double Volume

Optimum Growth™ 5L (3L Media)

vs Nalgene Nunc® 2.8L (1.5L)


**600%
INCREASE
IN PROTEIN**



**Optimum Growth™ 5L
1531mg/L
2X Volume & 3X Protein**

**Fernbach 2.8L
553mg/L**

**400%
INCREASE
IN PROTEIN**




**Optimum Growth™ 5L
430mg/L
2X Volume & 2X Protein**

**Fernbach 2.8L
266mg/L**

Conclusion

With Optimum Growth
2X Volume & 2-3X Protein



+



<



Fernbach 2.8L Flask
(1.5L Working Volume)

Fernbach 2.8L Flask
(1.5L Working Volume)

Optimum Growth™ 5L Flask
(3L Working Volume)

Recommended Shaker Speeds Optimum Growth™ 5 Liter Flasks

120RPM | 80RPM (2" Orbit)

