









**Tubes compatibility overview for Eppendorf Storage Boxes**

Box type	Height (inch/mm)	Grid	Ø Hole*	Box order no.	Tubes type
	2"/52.8 mm	10 x 10	12.5 mm	0030 140.508	Cryogenic tubes with internal thread, 1-2 mL
	2"/52.8 mm	9 x 9	13.15 mm	0030 140.516	Screw cap (cryogenic) tubes, 1-2 mL
	3"/76.2 mm	9 x 9	13.15 mm	0030 140.540	Screw cap (cryogenic) tubes, 3 mL
	4"/101.6 mm	9 x 9	13.15 mm	0030 140.567	Screw cap (cryogenic) tubes, 4-5 mL

\* For more information on dimensions see technical drawings available on the Storage Boxes product site at [www.eppendorf.com](http://www.eppendorf.com)