

Corning Incorporated  
Life SciencesRegistered  
ISO 9001:2008

# Product Description

**Catalog Number:** 1622-2**Product Description:** PYREX® Bottle, Specific Gravity, Gay Lussac, 2 mL  
with 7/15 Glass Stopper**Component Materials:**

Type I Class A - Low-Expansion Borosilicate Glass, per ASTM E438 (Standard Specification for Glasses in Laboratory Apparatus)

**Attributes:**Capacity - 2 ±0.1 mL  
Product Dimensions - 17 mm OD x 31 mm High without Stopper**Autoclaving:**

This product will withstand repeated standard autoclave cycles of 121°C at 15 psi for 15 minutes.

**Heating / Cooling:**

This product is primarily designed for a specific use that does not involve direct heating. In normal use, temperatures will be from -50°C to +125°C. This is below the temperature at which thermal failure could be a problem.

**Vacuum / Pressure:**

This product is not designed for and is not rated for vacuum or pressure applications.

**Special Features:**

Unadjusted for calibration in the laboratory.

**Sterilization:**

Shipped Not Sterile

**TSE/BSE:**

Materials of animal origin are not used in the manufacture of this product.

**Lot Number Designation:**

Each production lot is uniquely identified with an 8 digit lot number as follows: First 3 digits - Julian date, start of manufacturing; Next 2 digits - Year of manufacture; Last 3 digits - Batch identification.

Rev No: 1

**000 «Диаэм»**

Москва

ул. Магаданская, д. 7, к. 3 ■ тел./факс: (495) 745-0508 ■ sales@dia-m.ru

**www.dia-m.ru****С.-Петербург**+7 (812) 372-6040  
spb@dia-m.ru**Новосибирск**+7 (383) 328-0048  
nsk@dia-m.ru**Воронеж**+7 (473) 232-4412  
vrn@dia-m.ru**Йошкар-Ола**+7 (927) 880-3676  
nba@dia-m.ru**Красноярск**+7 (923) 303-0152  
krsk@dia-m.ru**Казань**+7 (843) 210-2080  
kazan@dia-m.ru**Ростов-на-Дону**+7 (863) 303-5500  
rnd@dia-m.ru**Екатеринбург**+7 (912) 658-7606  
ekb@dia-m.ru**Кемерово**+7 (923) 158-6753  
kemerovo@dia-m.ru**Армения**+7 (094) 01-0173  
armenia@dia-m.ru