












# Thermo Scientific™ Samco™ Exact Volume Transfer Pipette

- ▶ Polyethylene, non-sterile
- ▶ Useful for small samples or dispensing reagents
- ▶ Ideal for home & POC (Point of Care) testing
- ▶ Useful in clinical laboratories
- ▶ OEM applications
- ▶ Custom design options available to meet your specifications

**Cat. No.** **Dispense volume (µl)** **Stem Diameter (mm)** **Length (cm)** *Exact volume pipette lengths indicated below are close to actual size*

783NL	20	2.0	5.28	
784NL	40	2.3	5.92	
930NL	50	2.5	5.92	
785NL	60	2.5	6.48	
786NL	80	2.5	7.49	
787NL	100	2.5	8.30	
788NL	150	3.0	7.95	
941NL	200	3.0	8.89	
191NL	250	3.0	10.19	
966NL	300	3.6	9.50	
967NL	400	3.8	10.59	

500 pipets per pack in all sizes. 5000 pipets per case in all sizes.

**Thermo**  
SCIENTIFIC

**000 «Диаэм»**

Москва  
ул. Магаданская, д. 7, к. 3 ■ тел./факс: (495) 745-0508 ■ sales@dia-m.ru

**www.dia-m.ru**

**Новосибирск**  
пр. Академика  
Лаврентьева, д. 6/1  
тел.  
(383) 328-0048  
nsk@dia-m.ru

**Казань**  
ул. Парижской  
Коммуны, д. 6  
тел.  
(843) 210-2080  
kazan@dia-m.ru

**С.-Петербург**  
ул. Профессора  
Попова, д. 23  
тел.  
(812) 372-6040  
spb@dia-m.ru

**Ростов-  
на-Дону**  
пер. Семашко, д. 114  
тел.  
(863) 303-5500  
rnd@dia-m.ru

**Пермь**  
Представитель  
тел.  
(342) 202-2239  
perm@dia-m.ru

**Воронеж**  
Представитель  
тел.  
(473) 232-4412  
voronezh@dia-m.ru

**Армения**  
Представитель  
тел.  
(094) 01-0173  
armenia@dia-m.ru

**Узбекистан**  
Представитель  
тел.  
(90) 354-8569  
uz@dia-m.ru