







Eppendorf HeatSealer S100 and S200 Plate and Adapter Types

S100 (170 °C)	S200 (125 °C – 200 °C)	Suitable Plate Types
Plate Adapter, high profile 5391 070.034 	Plate Adapter, high profile 5392 070.011 	96-well Microplates 384-well Microplates 384-well PCR Plates Deepwell Plates 384/200µL 96/500µL
96-well PCR Plate Adapter 5391 070.026 	96-well PCR Plate Adapter 5392 070.038 	96-well PCR Plates unskirted semi-skirted skirted
Plate Adapter, low profile 5391 070.018 	Plate Adapter, low profile 5392 070.020 	Deepwell Plates 96/1000µL 96/2000µL

